


Whether you are transferring a patient from a bed, chair or floor, the portable Voyager will always give you the same dependability and reliability.

Our wide range of slings will help you in all types of situations.




**PORTABLE**

FEATURES	
Weight - kg (pounds)	5.5 kg (12 lbs.)
Lifting capacity	200 kg (440 lbs.)
Audible low-battery indicator	YES
Soft start and stop	YES
Emergency lowering device	MANUAL
Emergency stopping device	YES
CSA 601.1, UL 2601.1, EN 60601.1	YES


**ERGOFIT SLINGS**

With the specially designed BHM Medical slings, your patients should never have to sit on their sling again. All slings are designed to be installed and removed quickly without lifting the patient.




**HAMMOCK SLING**

For most transfers. Head-rest incorporated. Suitable for 80% of all patients. Excellent for bathing.




**QUICK-FIT SLING**

For general transfers. Perfectly suited for stocky or obese patients with strong hips or large thighs.



**HYGIENIC SLING**

Ideal for transferring to the toilet. It also permits an easy way to change diapers.



**WALKING SLING**

Helps to hold patients in a standing position. During walking exercises, can provide complete or partial support.

**SLING SIZES AND PRODUCT NUMBERS**

SIZE	SMALL	MEDIUM	LARGE
Suggested weight:	20 - 46 kg (45 - 100 lbs.)	46 - 95 kg (100 - 210 lbs.)	95 - 250 kg (210 - 550 lbs.)
Suggested height:	120 - 150 cm (4' - 4' 11")	151 - 180 cm (5' - 5' 11")	181 cm + (6' +)
Hammock	THA-S	THA-M	THA-L
Hammock 6	THA6-S	THA6-M	THA6-L
Quick-fit	TIR-S	TIR-M	TIR-L
Hygienic	THY-S	THY-M	THY-L
Walking	TEM-S	TEM-M	TEM-L

For a demonstration or more information:  
1-800-868-0441

2001 Tanguay Street  
Magog, Quebec, Canada J1X 5Y5  
Tel.: (819) 868-0441  
Fax: (819) 868-2249  
Internet: www.bhm-medical.com  
E-mail: bhm@bhm-medical.com



**PORTABLE**



Printed in Canada  
B 8001


**PORTABLE**

*Voyager*


**Voyager Portable: A track lift system that gives you all the versatility and flexibility you will ever need.**




We can help you with any track configuration. Call us for a free estimate.




*Integrated controls allow you to operate the lift with or without a handset.*



*A 24-volt battery allows you to perform all your tasks worry-free. When necessary, it recharges in less than 2 hours.*



*The carry bar can accommodate most types of slings.*



*Our high-strength track system easily installs onto any ceiling - new or old - while being discreet and esthetically pleasing.*



*Need a lift in more than one room? The Voyager only weighs 12 lbs. (battery included) and is designed to be carried.*

*Wonder how you can easily attach the lift to the track? Our unique "Quick Release" (patent pending) feature allows you to disengage the strap and quickly install or remove the lift from the track.*



*What if installing a ceiling track is not an option? Our semi-permanent rack gives you the same security as a permanent installation.*



*Combined with the Quick Release function, the "arm extension" gives people of any height the same ability to attach the lift to the track. At the same time unnecessary strain is avoided for the caregiver.*



*Our walking sling with the portable Voyager is a perfect match for those who need to be able to walk.*